

Examination Bed



Material:

The frame by steel coated. The quality attains national standard.

Introduction: PU surface with sponge inside.

Specifications:

- ◆ Outside size: L1900*W650*H460–900mm
- ◆ Inside size of mattress: L1900*W650mm
- ◆ Frame: Steel coated
- ◆ Surface: Equipped with 5cm soft sponge mattress, PU cover
- ◆ Legs: Bottoms of the legs have rubber pad, anti-skidding
- ◆ Feature: Waterproof, easy clean, durable
- ◆ Safe working load: 180kgs

