

## □ WSX - G1

合金铝楼梯担架  
ALUMINUM ALLOY STAIR CHAIR  
WSX-G合金铝楼梯担架采用优质  
铝合金材料,结构轻巧灵活,可折叠  
式结构便于运输携带,用于医院救  
护患者上下楼道之用.



## □ WSX - G1 吸塑板

合金铝楼梯担架  
ALUMINUM ALLOY STAIR CHAIR  
\* Thermoplastic Sitting and Back Part



- This article is an ideal facility for emergency staff to transfer patient to go up and down stair.
- Aluminum frame light and easy to handle.
- Two stationary rear wheels enable the patient to be moved over most surfaces with minimal effort.
- 3 pieces of restraint straps with quick-release buckle(see accessories items page)
- Foldable and save the space in ambulance car.

WSX - G1	承重 Capacity (kg)	净重 Net Weight (kg)	产品尺寸 Product Dimension (mm)	折叠尺寸 Folding Dimension (mm)
	159	8.5	690 x 510 x 910	180 x 510 x 910