

STAIR STRETCHER

Detailed Product Description:

This stretcher is designed for hospital emergency, emergency center and fire emergency, it can automatically slide down the stairs.

One medical staff can operate the stretcher, but better with another man to help

The height can be adjust according to the height of the patient.

Fasten the belts when patient seated, ensure the safety.

It can be folded, convenient for carry and storage.

Technical Features:

- 1) Product Size (L×W×H): 165×48.5×120cm
- 2) Folding Size (L×W×H): 109×48.5×25cm
- 3) Packing Size (1pc/carton): 116×58×32cm
- 4) N.W.: 33kg G.W.: 44kg
- 5) Loading Bearing: less than 180kg

