

# **BINOCULAR MICROSCOPE (XSZ-107BN)**

## **MULTI-PURPOSE BIOLOGICAL MICROSCOPE**



TC/XSP-N107



TC/XSZ-107BN



TC/XSZ-207



TC/XSP-106

	Specification	Code			
		TC/XSZ-N107	TC/XSZ-107BN	TC/XSZ-207	TC/XSP-106
<b>Viewing Head</b>	Sliding Binocular(45°),Rrotalatable(360°) Compensation Free Binocular(30°),Rrotalatable(360°)	●	●	●	●
<b>Nosepiece</b>	Quadruple	●	●	●	●
<b>Eyepiece</b>	WF10X(WF16X For Choice) H5X,H10X,H12.5X(H16X For Choice)	●	●	●	●
<b>Achromatic Objective</b>	DIN 4X, 10X, 40x(S), 100x(S, Oil)				
<b>Stage</b>	Double Layer Mechanocal Stage, Moving Range 75*45mm	140*140mm	140*140mm	140*140mm	130*125mm
<b>Focusing</b>	Coaxial Coarse Fine Adjustment, Focusing Rang 30mm Focusing Interval 0.002mm	●	●	●	●
<b>Condenser</b>	Abbe Condenser NA 1.25 With Iris Diaphragm & Filter	●	●	●	●
<b>Illumination Condenser</b>	Plain-concave Mirror Bult-in Alumination, Halogen Lamp 220v/110v, 6v, 20w	●	●	●	●