

EXAMINATION HALOGEN LIGHT

Halogen / Normal Examination Lights



HL-1220-PP



HL-1220-WM

Features:-

- * Non-magnetic Stainless Steel central tubes.
- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 700 mm.
- * Step down transformer 230V, 50W.
- * Halogen Bulb 12 V, 35 W, 10° angle of emission with UV block
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction.

Cat. No.	Description
HL-1220-PP	Halogen Examination Light, Adjustable Height, Mobile.
HL-1210-FH	Halogen Examination Light, Fixed Height, Mobile.
HL-1220-WM	Halogen Examination Light, Wall Mounted.



HL-1210-FH



EL-1230-PS

Features:-

- * Non-magnetic Stainless Steel central tubes.
- * Reinforced nylon base with 5 feet.
- * Mobile on twin wheel swivel castors Ø 50 mm .
- * Flexible goose neck arm 500 mm.
- * Normal Bulb 220 V, 60 W.
- * Input voltage : 220/230V 50-60Hz.
- * Knock down construction.

Cat. No.	Description
EL-1230-PS	Examination Light Fixed Height, Mobile.