



## **Emergency Trolley.**

- **Model No.: ET-63072B**
- **Specification (mm): 630 \* 470 \* 920**
- **Accessories & Details:**
  - **One-piece ABS plastic top board with raised-edge design, covered transparent soft plastic glass**
  - **Stainless steel guard rail, defibrillator shelf, sliding side shelf, storage box IV pole, dust basket**
  - **CPR board, British power outlet & hooks oxygen tank holder**
  - **Four Luxurious noiseless casters (two with brakes)**
  - **Centralized lock; 2 small drawers, 2 middle drawers and 1 big drawer (each inner with dividers can be organized easily and freely)**