



## **HOSPITAL BED 2 CRANK MANUAL**

### **Specifications:**

- **Outside size: L2150 \* W950 \* H500 mm**
- **Inside size of mattress base: L1950 \* W830 mm**
- **Head & Foot board: Detachable ABS Head & Foot board**
- **Guardrails: Collapsible aluminum guardrails**
- **Crank: 2 PCS Manual crank**
- **Wheels: 5" Castor with brakes**
- **Urine hooks: Urine hooks by steel chrome**
- **Safe working load: The quality attains national standard. It is equipped with silent castor, which is quiet and noiseless and able to support over 200 kgs.**