

PULSE OXIMETER FINGER TIP

H10 Finger Oximeter

H10 Finger Oximeter delivers a fast, reliable oximetry in an extremely handy, pocket-size solution by combining the monitor and the sensor into one unit. H10 is small in size; however it provides SpO₂, pulse rate measurement, bar graph and plethysmogram display.



Features:

- Small & lightweight design
- Dual color OLED display
- SpO₂, PR, Bar graph and Plethysmogram display
- 4 direction display with 6 display modes
- Easy-to-use
- 10-level adjustable brightness
- Low power consumption, automatic power off design
- Low-battery indicator
- Standard 2 AAA batteries supply



Single button operation

- Power on
- Pulse tone
- Display mode
- Screen brightness

